

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY A2 DPS

DATE 02/05/12 TIME 0800

INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D.

- HB* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution NA AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- HB* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- HB* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- HB* 4. Air blank completed.
- HB* 5. Calibration check completed. Test results 0.099 AC.
- HB* 6. Air blank completed.
- HB* 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III #638

INTOXILYZER 8000

Location: WICRENBURG PD

Serial Number: 89-000929

Core Version: 8105.49

02/05/2012

08:06:50

Standard Lot#: 25310100A1

Last Changed By: R. BYRNE, III #6386

DOS: R. BYRNE III #6386

AZ DPS

Operator: R. BYRNE III #6386

AZ DPS

Test	g/210L	Time
Air Blank	0.000	08:07:29
0.100 Cal Check	0.099	08:07:39
Air Blank	0.000	08:08:19

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 02/05/12 TIME 0811

INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG PD

AB 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

AB 1. Clock time check.

AB 2. Date check.

OPERATIONAL TESTS

AB 1. Alcohol-free subject test result 0.000 AC.

AB 2. Error recognition logic system functioning.

Not a Successfully Completed Test Sequence printed

AB 3. Proper sample recognition system.

Not a Successfully Completed Test Sequence printed

Deficient sample printed.

AB 4. Standard Calibration Check standard 0.100 AC. Result 0.098 AC.

Instrument is operating properly and accurately. YES AB NO

COMMENTS

LOT # 25310100A Exp 10/01/12

SIGNATURE Richard J. Byrne III *6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

02/05/2012

08:11:40

Standard Lot#: 25310100A1

Last Changed By: R. BYRNE, III #6386

OAS: R. BYRNE III #6386

AZ DPS

Operator: R. BYRNE III #6386

AZ DPS

Subject: SOAP,TEST

DOB: 08/08/1968

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:12:40
Diagnostic Test	Pass	08:13:15
Air Blank	0.000	08:13:42
0.100 Cal Check	0.099	08:14:03
Air Blank	0.000	08:14:33
Subject Test	0.000	08:14:50
Air Blank	ABT*	08:14:57
Air Blank	0.000	08:15:25

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

02/05/2012

08:15:35

Standard Lot#: 25310100A1

Last Changed By: R. BYRNE, III #6386

OAS: R. BYRNE III #6386

AZ DPS

Operator: R. BYRNE III #6386

AZ DPS

Subject: SOAP,TEST

DOB: 08/08/1968

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	08:16:11
Air Blank	RFI**	08:16:19

*Improper Sample

**RFI Detect

Not a Successfully Completed Test Sequence

INSTRUMENT# 8000
LOCATION: WILCHENBURG PD
Serial Number: 80-000929
Core Version: 8105.48
02/05/2012

08:16:34

Standard Lot#: 25310100A1
Last Changed By: R. BYRNE, JII #5386

QAS: R. BYRNE JII #6386
AZ DPS

Operator: R. BYRNE JII #6386
AZ DPS

Subject: SOAP TEST
DOB: 08/08/1968
Sex: M Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:17:33
Diagnostic Test	Pass	08:18:16
Air Blank	0.000	08:18:37
0.100 Cal Check	0.099	08:18:50
Air Blank	0.000	08:19:27
Subject Test	0.218*	08:22:05
Air Blank	0.000	08:23:07
Air Blank	0.000	08:23:34
Subject Test	0.204	08:23:51
Air Blank	IPG**	08:24:10
Air Blank	0.000	08:24:30

*Deficient Sample
**Improper Sample

REGISTERED IMPROPER
SAMPLE DURING AIR
BLANK

Was a Successfully Completed Test Sequence

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

02/05/2012

08:24:34

Standard Lot#: 25310100A1

Last Changed By: R. BYRNE, III

#6386

OAS: R. BYRNE III

#6386

AZ DPS

Operator: R. BYRNE III

#6386

AZ DPS

Subject: SQAP, TEST

DOB: 08/08/1968

Sex: M

Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	08:25:31
Diagnostic Test	Pass	08:26:06
Air Blank	0.000	08:26:34
0.100 Cal Check	0.098	08:26:55
Air Blank	0.000	08:27:24
Subject Test	0.257*	08:30:31
Air Blank	0.000	08:31:04
Air Blank	0.000	08:31:32
Subject Test	XXX**	08:31:52
Air Blank	0.000	08:32:28

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Richard J. Byrne III AGENCY AZ DPS
DATE 02/05/12 TIME 0839
INTOXILYZER SERIAL # 80-000929 DPS LOCATION Wickenburg P.D.

- RB* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- RB* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- RB* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- RB* 4. Air blank completed.
- RB* 5. Calibration check completed. Test results 0.098 AC.
- RB* 6. Air blank completed.
- RB* 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III #6386

HP0011222P 3000
Location: WICHITA/KR VP
Serial Number: 80-000029
Core Version: 9305108

02/05/2012 09:30:32

Standard Bolt: 25310100A1
Last Changed by: P. RYNER, III #6336

CAS: R. DYER #6336
32 DEC 2011

Operator: R. RYNER III #6336
AZ DES

Test	2/10/12	Time
Air Blank	0.000	06:39:00
0.100 Cal Check	0.098	08:24:21
Air Blank	0.000	08:24:53

Preventative Maintenance
Performed -- Peter Resce