

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 05/03/2012 TIME 0806

INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D.

- RB 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- RB 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- RB 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- RB 4. Air blank completed.
- RB 5. Calibration check completed. Test results 0.098 AC.
- RB 6. Air blank completed.
- RB 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard Byrne III #6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

05/03/2012

08:06:28

-----  
Standard Lot#: 25310100A1

Last Changed By: R. BYRNE, III #6386

-----  
QAS: R. BYRNE III #6386

AZ DPS

-----  
Operator: R. BYRNE III #6386

AZ DPS

-----  
Test                      g/210L                      Time  
-----  
Air Blank                      0.000                      08:07:01  
0.100 Cal Check                0.098                      08:07:22  
Air Blank                      0.000                      08:07:51

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS  
DATE 05/03/2012 TIME 0809  
INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D.

AS 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- AS 1. Clock time check.  
AS 2. Date check.

OPERATIONAL TESTS

- AS 1. Alcohol-free subject test result 0. 000 AC.  
AS 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
AS 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
AS 4. Standard Calibration Check standard 0. 120 AC. Result 0. 097 AC.

Instrument is operating properly and accurately. YES AS NO \_\_\_\_\_

COMMENTS

LOT # 25310100A1 Exp 10/01/2012

SIGNATURE Richard J. Byrne III #6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

05/03/2012

08:09:12

---

Standard Lot#: 25310100A1

Last Changed By: R. BYRNE, III #6386

---

QAS: R. BYRNE III #6386

AZ DPS

---

Operator: R. BYRNE III #6386

AZ DPS

---

Subject: SQAP, TEST

DOB: 08/08/1968

Sex: M

Weight: 185

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	08:10:10
Diagnostic Test	Pass	08:10:45
Air Blank	0.000	08:11:13
0.100 Cal Check	0.098	08:11:33
Air Blank	0.000	08:12:03
Subject Test	0.000	08:12:19
Air Blank	ABT*	08:12:26
Air Blank	0.000	08:12:54

\*Sequence Aborted

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: WICKENBURG PD  
Serial Number: 80-000929  
Core Version: 8105.48  
05/03/2012

08:13:13

---

Standard Lot#: 25310100A1  
Last Changed By: R. BYRNE, III #6386

---

OAS: R. BYRNE III #6386  
AZ DPS

---

Operator: R. BYRNE III #6386  
AZ DPS

---

Subject: SQAP,TEST  
DOB: 08/08/1968  
Sex: M Weight: 185

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	IPS*	08:13:49
Air Blank	RFI**	08:14:01

\*Improper Sample  
\*\*RFI Detect

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

05/03/2012

08:14:12

---

Standard Lot#: 25310100A1

Last Changed By: R. BYRNE, III #6386

---

QAS: R. BYRNE III #6386

AZ DPS

---

Operator: R. BYRNE III #6386

AZ DPS

---

Subject: SQAP, TEST

DOB: 08/08/1968

Sex: M

Weight: 185

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	08:15:05
Diagnostic Test	Pass	08:15:40
Air Blank	0.000	08:16:08
0.100 Cal Check	0.098	08:16:29
Air Blank	0.000	08:16:58
Subject Test	0.268*	08:20:06
Air Blank	0.000	08:20:38
Air Blank	0.000	08:21:06
Subject Test	XXX**	08:21:22
Air Blank	0.000	08:21:56

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 05/03/2012 TIME 0823

INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D.

- (AB) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution, 100 AC.  
OR  
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (AB) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (AB) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (AB) 4. Air blank completed.
- (AB) 5. Calibration check completed. Test results 0.097 AC.
- (AB) 6. Air blank completed.
- (AB) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III #6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

05/03/2012

08:23:13

-----  
Standard Lot#: 25310100A1

Last Changed By: R. BYRNE, III #6386

-----  
QAS: R. BYRNE III #6386

AZ DPS

-----  
Operator: R. BYRNE III #6386

AZ DPS

-----  
Test                      g/210L                      Time  
-----  
Air Blank                      0.000                      08:23:44  
0.100 Cal Check                0.097                      08:24:05  
Air Blank                      0.000                      08:24:34

Preventative Maintenance  
Performed -- Timer Reset