

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE AGENCY AZ DPS

DATE 10/16/2013 TIME 1056

INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D.

- (FB) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
Pour a standard alcohol concentration solution OR NA AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (FB) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (FB) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (FB) 4. Air blank completed.
- (FB) 5. Calibration check completed. Test results 0.098 AC.
- (FB) 6. Air blank completed.
- (FB) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III #6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

10/16/2013

10:56:42

-----  
Standard Lot#: 15712100A1

Last Changed By: R. BYRNE III #6386  
-----

QAS: R. BYRNE III #6386

AZ DPS  
-----

Operator: R. BYRNE III #6386

AZ DPS  
-----

Test	g/210L	Time
Air Blank	0.000	10:57:12
0.100 Cal Check	0.098	10:57:33
Air Blank	0.000	10:58:03

Preventative Maintenance  
Performed -- Timer Reset

EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 10/16/2013 TIME 1100

INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D.

AB 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

AB 1. Clock time check.

AB 2. Date check.

OPERATIONAL TESTS

AB 1. Alcohol-free subject test result 0. 000 AC.

AB 2. Error recognition logic system functioning.

Not a Successfully Completed Test Sequence printed

AB 3. Proper sample recognition system.

Not a Successfully Completed Test Sequence printed

Deficient sample printed.

AB 4. Standard Calibration Check standard 0. 100 AC. Result 0. 078 AC.

Instrument is operating properly and accurately. YES AB NO \_\_\_\_\_

COMMENTS

LOT # 15712100A1 Exp 08/01/2014

SIGNATURE Richard J. Byrne III #6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

10/16/2013

11:00:25

---

Standard Lot#: 15712100A1

Last Changed By: R. BYRNE III #6386

---

QAS: R. BYRNE III #6386

AZ DPS

---

Operator: R. BYRNE III #6386

AZ DPS

---

Subject: SQAP,TEST

DOB: 08/08/1968

Sex: M

Weight: 185

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	11:01:32
Diagnostic Test	Pass	11:02:07
Air Blank	0.000	11:02:35
0.100 Cal Check	0.099	11:02:56
Air Blank	0.000	11:03:25
Subject Test	0.000	11:03:53
Air Blank	0.000	11:04:20
Wait	ABT*	
Air Blank	0.000	11:04:56

\*Sequence Aborted

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000  
Location: WICKENBURG PD  
Serial Number: 80-000929  
Core Version: 8105.48  
10/16/2013

11:05:05

Standard Lot#: 15712100A1  
Last Changed By: R. BYRNE III #6386

QAS: R. BYRNE III #6386  
AZ DPS

Operator: R. BYRNE III #6386  
AZ DPS

Subject: SQAP, TEST  
DOB: 08/08/1968  
Sex: M Weight: 185

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	11:05:42
Air Blank	RFI**	11:05:50

\*Improper Sample

\*\*RFI Detect

**Not a Successfully Completed Test Sequence**

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

10/16/2013

11:05:59

---

Standard Lot#: 15712100A1

Last Changed By: R. BYRNE III #6386

---

QAS: R. BYRNE III #6386

AZ DPS

---

Operator: R. BYRNE III #6386

AZ DPS

---

Subject: SQAP,TEST

DOB: 08/08/1968

Sex: M Weight: 185

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	11:06:57
Diagnostic Test	Pass	11:07:32
Air Blank	0.000	11:07:59
0.100 Cal Check	0.099	11:08:20
Air Blank	0.000	11:08:49
Subject Test	0.094*	11:11:57
Air Blank	0.000	11:12:29
Air Blank	0.000	11:12:56
Subject Test	XXX**	11:13:20
Air Blank	0.000	11:13:54

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

**Not a Successfully Completed Test Sequence**

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)  
ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS  
DATE 10/16/2013 TIME 1114  
INTOXILYZER SERIAL # 80-000929 LOCATION WICKENBURG P.D.

- AB 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.  
OR  
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- BB 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- BB 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- BB 4. Air blank completed.
- BB 5. Calibration check completed. Test results 0. 098 AC.
- BB 6. Air blank completed.
- BB 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III # 6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

10/16/2013

11:14:11

-----  
Standard Lot#: 15712100A1

Last Changed By: R. BYRNE III #6386  
-----

QAS: R. BYRNE III #6386

AZ DPS  
-----

Operator: R. BYRNE III #6386

AZ DPS  
-----

Test	g/210L	Time
Air Blank	0.000	11:14:43
0.100 Cal Check	0.098	11:15:04
Air Blank	0.000	11:15:33

Preventative Maintenance  
Performed -- Timer Reset