

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 07/12/2014 TIME 0958

INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D.

- RB* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution 100 AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- RB* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- RB* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- RB* 4. Air blank completed.
- RB* 5. Calibration check completed. Test results 0.028 AC.
- RB* 6. Air blank completed.
- RB* 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE *Richard J. Byrne III* #6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

07/12/2014

09:58:41

Standard Lot#: 15712100A1

Last Changed By: R. BYRNE III #6386

QAS: R. BYRNE III #6386

AZ DPS

Operator: R. BYRNE III #6386

AZ DPS

Test g/210L Time

Air Blank 0.000 09:59:23
0.100 Cal Check 0.098 09:59:44
Air Blank 0.000 10:00:13

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 02/12/2014 TIME 1010

INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D.

- RB* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 1100 AC.
OR
Pour a standard alcohol concentration solution 11/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- RB* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- RB* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- RB* 4. Air blank completed.
- RB* 5. Calibration check completed. Test results 0.095 AC.
- RB* 6. Air blank completed.
- RB* 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE *Richard J. Byrne III* # 6386

INTOXILYZER 8000

Location: WICKENBURG PD

Serial Number: 80-000929

Core Version: 8105.48

07/12/2014

10:10:12

Standard Lot#: 15712100A1

Last Changed By: R. BYRNE III #6386

QAS: R. BYRNE III #6386

AZ DPS

Operator: R. BYRNE III #6386

AZ DPS

Test g/210L Time

Air Blank 0.000 10:10:44
0.100 Cal Check 0.097 10:11:05
Air Blank 0.000 10:11:34

Preventative Maintenance
Performed -- Timer Reset

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 07/21/2014 TIME 1012

INTOXILYZER SERIAL # 80-000929 DPS LOCATION WICKENBURG P.D.

AB 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- AB 1. Clock time check.
AB 2. Date check.

OPERATIONAL TESTS

- AB 1. Alcohol-free subject test result 0.00 AC.
AB 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
AB 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
AB 4. Standard Calibration Check standard 0.600 AC. Result 0.098 AC.

Instrument is operating properly and accurately. YES AB NO _____

COMMENTS LOT # 15712100A1 EXP 08/01/2014

ANNUAL RECERT

SIGNATURE Richard J. Byrne III #6386

INTOXILYZER 8000
Location: WICKENBURG PD
Serial Number: 80-000929
Core Version: 8105.48
07/12/2014 10:12:45

Standard Lot#: 15712100A1
Last Changed By: R. BYRNE III #6386

QAS: R. BYRNE III #6386
AZ DPS

Operator: R. BYRNE III #6386
AZ DPS

Subject: SQAP,TEST
DOB: 08/08/1968
Sex: M Weight: 185

15 Minute Deprivation Period? Yes

| Test | g/210L | Time |
|-----------------|--------|----------|
| Air Blank | 0.000 | 10:13:42 |
| Diagnostic Test | Pass | 10:14:17 |
| Air Blank | 0.000 | 10:14:44 |
| 0.100 Cal Check | 0.098 | 10:15:05 |
| Air Blank | 0.000 | 10:15:35 |
| Subject Test | 0.000 | 10:15:51 |
| Air Blank | 0.000 | 10:16:19 |
| Wait | ABT* | |
| Air Blank | 0.000 | 10:16:55 |

*Sequence Aborted

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: WICKENBURG PD
Serial Number: 80-000929
Core Version: 8105.48
07/12/2014 10:17:00

Standard Lot#: 15712100A1
Last Changed By: R. BYRNE III #6386

QAS: R. BYRNE III #6386
AZ DPS

Operator: R. BYRNE III #6386
AZ DPS

Subject: SQAP,TEST
DOB: 08/08/1968
Sex: M Weight: 185

15 Minute Deprivation Period? Yes

| Test | g/210L | Time |
|-----------|--------|----------|
| Air Blank | IPS* | 10:17:39 |
| Air Blank | RFI** | 10:17:46 |

*Improper Sample
**RFI Detect

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: WICKENBURG PD
Serial Number: 80-000929
Core Version: 8105.48
07/12/2014

10:18:05

Standard Lot#: 15712100A1
Last Changed By: R. BYRNE III #6386

QAS: R. BYRNE III #6386
AZ DPS

Operator: R. BYRNE III #6386
AZ DPS

Subject: SQAP,TEST
DOB: 08/08/1968
Sex: M Weight: 185

15 Minute Deprivation Period? Yes

| Test | g/210L | Time |
|-----------------|--------|----------|
| Air Blank | 0.000 | 10:19:05 |
| Diagnostic Test | Pass | 10:19:40 |
| Air Blank | 0.000 | 10:20:07 |
| 0.100 Cal Check | 0.098 | 10:20:28 |
| Air Blank | 0.000 | 10:20:58 |
| Subject Test | 0.206* | 10:24:05 |
| Air Blank | 0.000 | 10:24:39 |
| Air Blank | 0.000 | 10:25:06 |
| Subject Test | XXX** | 10:26:16 |
| Air Blank | 0.000 | 10:26:49 |

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST RICHARD J. BYRNE III AGENCY AZ DPS

DATE 07/14/2014 TIME 1030

INTOXILYZER SERIAL # 80-000929 LOCATION WICKENBURG P.D.

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 098 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Richard J. Byrne III # 63236

