

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST SCARIM, ROBERT G AGENCY COTTONWOOD PD
DATE 5/19/08 TIME 1935
INTOXILYZER SERIAL # 80-000933 LOCATION COTTONWOOD PD

- (-1) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (M) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (M) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (M) 4. Air blank completed.
- (M) 5. Calibration check completed. Test results 0.098 AC.
- (M) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE ^{OK} [Signature] 592

DPS Form Exh G-2 (Rev 05-01)

5/30/08 Copy to PDD (AK)

08-004846

INTOXILYZER 8000

Location: COTTONWOOD PD

Serial Number: 80-000933

05/19/2008

19:37:49

Standard Lot#: 7179011


Last Changed By: R. SCARIM #592

QAS: R. SCARIM #592

COTTONWOOD PD

Test	g/210L	Time
Air Blank	0.000	19:38:16
Cal Check	0.098	19:38:35
Air Blank	0.000	19:39:05

OFC



Operator's Signature

08-004846