

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)
ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Ofc M. WATSON AGENCY COTTONWOOD PD
DATE 09/18/2008 TIME 0042 HRS
INTOXILYZER SERIAL # 80-000933 LOCATION COTTONWOOD PD

- (MWS) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (MWS) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (MWS) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (MWS) 4. Air blank completed.
- (MWS) 5. Calibration check completed. Test results 0.098 AC.
- (MWS) 6. Air blank completed.
- (MWS) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE Ofc [Signature] #6109

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number: 80-000933
09/18/2008 00:42:09

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Test	g/210L	Time
Air Blank	0.000	00:42:37
Cal Check	0.098	00:42:55
Air Blank	0.000	00:43:25

M. Watson #609

Operator's Signature