

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST OFF M. WATSON AGENCY COTTONWOOD PD
DATE 04/25/2009 TIME 2001
INTOXILYZER SERIAL # 80-000933 LOCATION COTTONWOOD PD

- (MWS) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (MWS) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
(MWS) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
(MWS) 4. Air blank completed.
(MWS) 5. Calibration check completed. Test results 0.099 AC.
(MWS) 6. Air blank completed.
(MWS) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE OFF M. WATSON #609


INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number: 80-000933
04/25/2009

20:01:32

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Test	g/210L	Time
Air Blank	0.000	20:02:00
Cal Check	0.099	20:02:19
Air Blank	0.000	20:02:48

OFC  *A609*

Operator's Signature

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST OFC M. WATSON AGENCY COTTONWOOD PD
DATE 04/25/09 TIME 2003
INTOXILYZER SERIAL # 80-000933 LOCATION COTTONWOOD PD

~~()~~ 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- ~~()~~ 1. Clock time check.
~~()~~ 2. Date check.

OPERATIONAL TESTS

- ~~()~~ 1. Alcohol-free subject test result 0.000 AC.
~~()~~ 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
~~()~~ 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
() 4. Standard Calibration Check standard 0. _____ AC. Result 0. _____ AC.

Instrument is operating properly and accurately. YES _____ NO _____

COMMENTS

RAJ PROPER SAMPLE RECOGNITION TWO TIMES WITH RANGE EXCEEDED

SIGNATURE OFC  #609

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:03:36

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	20:04:40
Diagnostic Test	Pass	20:05:08
Air Blank	0.000	20:05:35
0.100 Cal Check	0.099	20:05:54
Air Blank	0.000	20:06:24
Subject Test	0.000	20:06:49
Air Blank	0.000	20:07:16
Five Minute Wait	ABT*	
Air Blank	0.000	20:07:58

*Sequence Aborted

Not a Successfully Completed Test Sequence



INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:03:36

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP, TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	20:04:40
Diagnostic Test	Pass	20:05:08
Air Blank	0.000	20:05:35
0.100 Cal Check	0.099	20:05:54
Air Blank	0.000	20:06:24
Subject Test	0.000	20:06:49
Air Blank	0.000	20:07:16
Five Minute Wait	ABT*	
Air Blank	0.000	20:07:58

*Sequence Aborted

Not a Successfully Completed Test Sequence



INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:03:36

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	20:04:40
Diagnostic Test	Pass	20:05:08
Air Blank	0.000	20:05:35
0.100 Cal Check	0.099	20:05:54
Air Blank	0.000	20:06:24
Subject Test	0.000	20:06:49
Air Blank	0.000	20:07:16
Five Minute Wait	ABT*	
Air Blank	0.000	20:07:58

*Sequence Aborted

Not a Successfully Completed Test Sequence



INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:08:17

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	20:08:53
Air Blank	0.000	20:09:21

*Improper Sample

Not a Successfully Completed Test Sequence



INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:08:17

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	20:08:53
Air Blank	0.000	20:09:21

*Improper Sample

Not a Successfully Completed Test Sequence



INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:08:17

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP, TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	20:08:53
Air Blank	0.000	20:09:21

*Improper Sample

Not a Successfully Completed Test Sequence



INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:10:01

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	20:10:35
Air Blank	0.000	20:11:04

*RFI Detect

Not a Successfully Completed Test Sequence



INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:10:01

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	20:10:35
Air Blank	0.000	20:11:04

*RFI Detect

Not a Successfully Completed Test Sequence



INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:10:01

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	20:10:35
Air Blank	0.000	20:11:04

*RFI Detect

Not a Successfully Completed Test Sequence



INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:11:19

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	20:12:09
Diagnostic Test	Pass	20:12:37
Air Blank	0.000	20:13:04
0.100 Cal Check	0.099	20:13:23
Air Blank	0.000	20:13:52
Subject Test	0.243*	20:17:01
Air Blank	0.000	20:17:33
Air Blank	0.000	20:18:00
Subject Test	RGE**	20:18:20
Air Blank	0.000	20:18:56

*Deficient Sample

**Range Exceeded

Not a Successfully Completed Test Sequence



INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:11:19

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	20:12:09
Diagnostic Test	Pass	20:12:37
Air Blank	0.000	20:13:04
0.100 Cal Check	0.099	20:13:23
Air Blank	0.000	20:13:52
Subject Test	0.243*	20:17:01
Air Blank	0.000	20:17:33
Air Blank	0.000	20:18:00
Subject Test	RGE**	20:18:20
Air Blank	0.000	20:18:56

*Deficient Sample

**Range Exceeded

Not a Successfully Completed Test Sequence



INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:11:19

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	20:12:09
Diagnostic Test	Pass	20:12:37
Air Blank	0.000	20:13:04
0.100 Cal Check	0.099	20:13:23
Air Blank	0.000	20:13:52
Subject Test	0.243*	20:17:01
Air Blank	0.000	20:17:33
Air Blank	0.000	20:18:00
Subject Test	RGE**	20:18:20
Air Blank	0.000	20:18:56

*Deficient Sample

**Range Exceeded

Not a Successfully Completed Test Sequence



INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:19:16

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	20:20:06
Diagnostic Test	Pass	20:20:33
Air Blank	0.000	20:21:01
0.100 Cal Check	0.098	20:21:19
Air Blank	0.000	20:21:49
Subject Test	RGE*	20:24:05
Air Blank	0.000	20:24:42

*Range Exceeded

Not a Successfully Completed Test Sequence



INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:19:16

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD


Subject: SQAP, TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	20:20:06
Diagnostic Test	Pass	20:20:33
Air Blank	0.000	20:21:01
0.100 Cal Check	0.098	20:21:19
Air Blank	0.000	20:21:49
Subject Test	RGE*	20:24:05
Air Blank	0.000	20:24:42

*Range Exceeded

Not a Successfully Completed Test Sequence



INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:19:16

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	20:20:06
Diagnostic Test	Pass	20:20:33
Air Blank	0.000	20:21:01
0.100 Cal Check	0.098	20:21:19
Air Blank	0.000	20:21:49
Subject Test	RGE*	20:24:05
Air Blank	0.000	20:24:42

*Range Exceeded

Not a Successfully Completed Test Sequence



INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:24:55

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	20:25:47
Diagnostic Test	Pass	20:26:15
Air Blank	0.000	20:26:42
0.100 Cal Check	0.097	20:27:01
Air Blank	0.000	20:27:30
Subject Test	0.286*	20:30:39
Air Blank	0.000	20:31:11
Air Blank	0.000	20:31:39
Subject Test	XXX**	20:32:11
Air Blank	0.000	20:32:45

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:24:55

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	20:25:47
Diagnostic Test	Pass	20:26:15
Air Blank	0.000	20:26:42
0.100 Cal Check	0.097	20:27:01
Air Blank	0.000	20:27:30
Subject Test	0.286*	20:30:39
Air Blank	0.000	20:31:11
Air Blank	0.000	20:31:39
Subject Test	XXX**	20:32:11
Air Blank	0.000	20:32:45

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
04/25/2009

20:24:55

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1901
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	20:25:47
Diagnostic Test	Pass	20:26:15
Air Blank	0.000	20:26:42
0.100 Cal Check	0.097	20:27:01
Air Blank	0.000	20:27:30
Subject Test	0.286*	20:30:39
Air Blank	0.000	20:31:11
Air Blank	0.000	20:31:39
Subject Test	XXX**	20:32:11
Air Blank	0.000	20:32:45

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST OFF M. WATSON AGENCY COTTONWOOD PD
DATE 04/25/2009 TIME 2033
INTOXILYZER SERIAL # 80-000933 LOCATION COTTONWOOD PD

1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
4. Air blank completed.
5. Calibration check completed. Test results 0. 086 AC.
6. Air blank completed.
7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE OFF [Signature] #609


INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number: 80-000933
04/25/2009

20:33:47

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Test	g/210L	Time
Air Blank	0.000	20:34:15
Cal Check	0.098	20:34:34
Air Blank	0.000	20:35:03

 #609

Operator's Signature