

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. WATSON AGENCY COTTONWOOD PD
DATE 07/21/09 TIME 0958
INTOXILYZER SERIAL # 80-000933 LOCATION COTTONWOOD PD

- W/W* 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- W/W* 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- W/W* 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- W/W* 4. Air blank completed.
- W/W* 5. Calibration check completed. Test results 0.098 AC.
- W/W* 6. Air blank completed.
- W/W* 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE OKC [Signature] #609

INTOXILYZER 8000

Location: COTTONWOOD PD

Serial Number: 80-000933

07/21/2009

09:58:32

Standard Lot#: 804301I

Last Changed By: R. SCARIM #592

QAS: M. WATSON #609

COTTONWOOD PD

Test	g/210L	Time
Air Blank	0.000	09:58:59
Cal Check	0.098	09:59:18
Air Blank	0.000	09:59:48

MFC [Signature] #609

Operator's Signature

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST M. WATSON AGENCY COTTONWOOD PD
DATE 07/21/09 TIME 1001
INTOXILYZER SERIAL # 80-000933 LOCATION COTTONWOOD PD

~~FAIL~~ 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- ~~FAIL~~ 1. Clock time check.
~~FAIL~~ 2. Date check.

OPERATIONAL TESTS

- ~~FAIL~~ 1. Alcohol-free subject test result 0. 200 AC.
~~FAIL~~ 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
~~FAIL~~ 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
~~FAIL~~ 4. Standard Calibration Check standard 0. 100 AC. Result 0. 098 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

SIGNATURE OFC [Signature] #609

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
07/21/2009

10:01:04

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD


Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1900
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	10:02:02
Diagnostic Test	Pass	10:02:29
Air Blank	0.000	10:02:57
0.100 Cal Check	0.099	10:03:15
Air Blank	0.000	10:03:45
Subject Test	0.000	10:04:05
Air Blank	0.000	10:04:33
Five Minute Wait	ABT*	
Air Blank	0.000	10:05:13

*Sequence Aborted



Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
07/21/2009

10:01:04

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD


Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1900
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	10:02:02
Diagnostic Test	Pass	10:02:29
Air Blank	0.000	10:02:57
0.100 Cal Check	0.099	10:03:15
Air Blank	0.000	10:03:45
Subject Test	0.000	10:04:05
Air Blank	0.000	10:04:33
Five Minute Wait	ABT*	
Air Blank	0.000	10:05:13

*Sequence Aborted



Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
07/21/2009

10:01:04

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1900
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	10:02:02
Diagnostic Test	Pass	10:02:29
Air Blank	0.000	10:02:57
0.100 Cal Check	0.099	10:03:15
Air Blank	0.000	10:03:45
Subject Test	0.000	10:04:05
Air Blank	0.000	10:04:33
Five Minute Wait	ABT*	
Air Blank	0.000	10:05:13

*Sequence Aborted



Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
07/21/2009

10:06:07

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1900
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	10:06:51
Air Blank	0.000	10:07:19

*Improper Sample



Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
07/21/2009

10:06:07

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP, TEST
DOB: 01/01/1900
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	10:06:51
Air Blank	0.000	10:07:19

*Improper Sample



Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
07/21/2009

10:06:07

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1900
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	10:06:51
Air Blank	0.000	10:07:19

*Improper Sample



Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
07/21/2009

10:08:17

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1900
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	10:09:02
Air Blank	0.000	10:09:30

*RFI Detect



Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
07/21/2009

10:08:17

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1900
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	10:09:02
Air Blank	0.000	10:09:30

*RFI Detect



Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
07/21/2009

10:08:17

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1900
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	10:09:02
Air Blank	0.000	10:09:30

*RFI Detect



Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
07/21/2009

10:10:09

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1900
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	10:10:59
Diagnostic Test	Pass	10:11:26
Air Blank	0.000	10:11:53
0.100 Cal Check	0.098	10:12:12
Air Blank	0.000	10:12:42
Subject Test	0.015*	10:15:50
Air Blank	0.000	10:16:22
Air Blank	0.000	10:16:49
Subject Test	XXX**	10:17:22
Air Blank	0.000	10:17:55

*Deficient Sample

**Invalid Sample

Begin new deprivation period



Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
07/21/2009

10:10:09

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1900
Sex: M Weight: 111


15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	10:10:59
Diagnostic Test	Pass	10:11:26
Air Blank	0.000	10:11:53
0.100 Cal Check	0.098	10:12:12
Air Blank	0.000	10:12:42
Subject Test	0.015*	10:15:50
Air Blank	0.000	10:16:22
Air Blank	0.000	10:16:49
Subject Test	XXX**	10:17:22
Air Blank	0.000	10:17:55

*Deficient Sample

**Invalid Sample

Begin new deprivation period



Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number : 80-000933
07/21/2009 10:10:09

Standard Lot#: 804301I
Last Changed By: R. SCARIM #592

QAS: M. WATSON #609
COTTONWOOD PD

Operator: M. WATSON #609
COTTONWOOD PD

Subject: SQAP,TEST
DOB: 01/01/1900
Sex: M Weight: 111

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	10:10:59
Diagnostic Test	Pass	10:11:26
Air Blank	0.000	10:11:53
0.100 Cal Check	0.098	10:12:12
Air Blank	0.000	10:12:42
Subject Test	0.015*	10:15:50
Air Blank	0.000	10:16:22
Air Blank	0.000	10:16:49
Subject Test	XXX**	10:17:22
Air Blank	0.000	10:17:55

*Deficient Sample

**Invalid Sample

Begin new deprivation period



Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST M. WATSON AGENCY COTTON WOOD PD
DATE 07/21/09 TIME 1019
INTOXILYZER SERIAL # 80-000933 LOCATION COTTON WOOD PD

- (S/W) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (S/W) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
(S/W) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
(S/W) 4. Air blank completed.
(S/W) 5. Calibration check completed. Test results 0. 098 AC.
(S/W) 6. Air blank completed.
(S/W) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE DFC [Signature] #609

INTOXILYZER 8000

Location: COTTONWOOD PD

Serial Number: 80-000933

07/21/2009

10:19:09

Standard Lot#: 804301I


Last Changed By: R. SCARIM #592

QAS: M. WATSON #609

COTTONWOOD PD

Test g/210L Time

Air Blank 0.000 10:19:37
Cal Check 0.098 10:19:55
Air Blank 0.000 10:20:25



Operator's Signature

PROPERTY INVOICE AND RECEIPT

NOTE: FILL IN EVERY SPACE UNLESS THE INFORMATION IS NOT AVAILABLE PRINT LEGIBLY IN INK.		TYPE OF PROPERTY (Circle One Only)		DATE AND TIME PROPERTY SEIZED	
		<u>EVIDENCE</u>	SAFEKEEPING	FOUND	
NAME OF VICTIM OR OWNER (LAST, FIRST, MIDDLE) <i>COTTONWOOD PD</i>		ADDRESS		TELEPHONE	D.R. NUMBER <i>1020 H25</i>
NAME OF PERSON ARRESTED, OR SUSPECT (LAST, FIRST, MIDDLE)		ADDRESS		TELEPHONE	D.O.B. <i>609-0016773</i>
NAME OF PERSON ARRESTED, OR SUSPECT (LAST, FIRST, MIDDLE)		ADDRESS		TELEPHONE	D.O.B.
NAME OF FINDER OR PERSON REPORTING FOUND PROPERTY		ADDRESS		TELEPHONE	CITATION #
CHARGES				FINDER WILL CLAIM PROPERTY	
				YES	<u>NO</u>
				YES	<u>NO</u>

NOTIFICATION FOR RETURN OF PERSONAL PROPERTY

You are responsible for claiming your personal property within 90 days from the above date. The property will be disposed of by the Cottonwood Police Department, if no claim is made prior to that date.

ADVANCE AUTHORIZATION IS REQUIRED IF YOU WISH TO DESIGNATE ANOTHER PERSON TO RECEIVE YOUR PROPERTY.

SIGNATURE

BRIEF SYNOPSIS OF OCCURRENCE: (Also use this space for additional suspects or victims)
30/90 DAY QAS

ITEMIZE AND DESCRIBE PROPERTY

LIST ONLY ONE ITEM PER LINE. NUMBER EACH ITEM ON THE FORM. - Where possible, place each item in plastic bag, and attach property tag to outside of bag.

ITEM	DESCRIPTION	SERIAL #	PROPERTY #
1	<i>30/90 DAY QAS ENTOX LYZER 3000</i>		<i>70449</i>

Placed in Property Locker # *3* Date *7/21/09* Time *1020*

OFFICER: *M. WATSON* SERIAL: *609* PAGE *1* OF *1*

DISPOSITION OF PROPERTY:
Property Items Were:

1. Returned to Owner	3. Sold at Auction	5. Left at Court
2. Destroyed	4. T.O.T. City	6. Unable to Locate

Received By <i>[Signature]</i>	Date <i>7/21/09</i>	Time <i>1430</i>	Remarks
Property Location <i>DWF</i>			