

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST SCARIM, ROGER G. AGENCY COTTONWOOD PD  
DATE 7/30/09 TIME 1646  
INTOXILYZER SERIAL # 80-000933 LOCATION COTTONWOOD PD

- (M) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (M) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (M) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (M) 4. Air blank completed.
- (M) 5. Calibration check completed. Test results 0.097 AC.
- (M) 6. Air blank completed.
- (M) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE of [Signature] 592

INTOXILYZER 8000  
Location: COTTONWOOD PD  
Serial Number: 80-000933  
07/30/2009

16:47:55

-----  
Standard Lot#: 585219  
Last Changed By: R. SCARIM #592  
-----

QAS: R. SCARIM #592  
COTTONWOOD PD  
-----

Test	g/210L	Time
Air Blank	0.000	16:48:22
Cal Check	0.097	16:48:41
Air Blank	0.000	16:49:11

QFC



Operator's Signature

\* BOTTLE # 585219 INSTALLED

39-007049

