



EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST SCARON, ROGER G. AGENCY COTTONWOOD PD  
DATE 2/27/09 TIME 1700  
INTOXILYZER SERIAL # 30-000933 LOCATION COTTONWOOD PD

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0.100 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

INTOXILYZER 8000

Location: COTTONWOOD PD

Serial Number: 80-000933

02/27/2009

17:00:50


-----  
Standard Lot#: 804301I

Last Changed By: R. SCARIM #592  
-----

QAS: R. SCARIM #592

COTTONWOOD PD  
-----

Test	g/210L	Time
Air Blank	0.000	17:01:17
Cal Check	0.100	17:01:36
Air Blank	0.000	17:02:05

*ofc*  
  
592

Operator's Signature