

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST SCARIM, ROGER G. AGENCY COTTONWOOD PD  
DATE 4/4/18 TIME 1329  
INTOXILYZER SERIAL # 80-000933 LOCATION COTTONWOOD PD

- (M) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.  
OR  
Pour a standard alcohol concentration solution \_\_\_\_\_ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (M) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"  
(M) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.  
(M) 4. Air blank completed.  
(M) 5. Calibration check completed. Test results 0.100 AC.  
(M) 6. Air blank completed.  
(M) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE OR [Signature] 592

DPS Form Exh G-2 (Rev 05-01)

BOTTLE LOT # 0261810091 ADDED

INTOXILYZER 8000  
Instrument Initialization  
13:01:05 04/04/2018

INTOXILYZER 8000  
Location: COTTONWOOD PD  
Serial Number: 80-000933  
Core Version: 8105.48  
04/04/2018 13:31:26

Standard Lot#: 0261810091  
Last Changed By: R. SCARIM #592

GAS: R. SCARIM #592  
COTTONWOOD PD

Operator: R. SCARIM #592  
COTTONWOOD PD

Test	g/210L	Time
Air Blank	0.000	13:31:59
0.100 Cal Check	0.100	13:32:20
Air Blank	0.000	13:32:49

Preventative Maintenance  
Performed -- Timer Reset