

DIR 17-011168

COPY

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST S. Blakley Sr. AGENCY Yavapai County SO

DATE 3-30-17 TIME 21:16:08

INTOXILYZER SERIAL # 80-001019 LOCATION Prescott YCSO

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution NA AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0. 098 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

DPS Form Exh G-2 (Rev 05-01)

Tank Serial # 27016100A1

Exp 11-5-2018

Press 900 PSI

RECEIVED

APR 04 2017

YAVAPAI COUNTY ATTORNEY

DR 17-011168

INTOXILYZER 8000

Location: YCSO PRESCOTT

Serial Number: 80-001019

Core Version: 8105.48

03/30/2017

21:16:08

Standard Lot#: 27016100A1

Last Changed By: S. BLAKLEY SR #L6861

QAS: S. BLAKLEY SR #L6861

YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	21:16:42
0.100 Cal Check	0.098	21:17:03
Air Blank	0.000	21:17:32

Preventative Maintenance
Performed -- Timer Reset

DR 17-011168

COPY

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST S. Blankley Sr. AGENCY Yavapai County SO

DATE 3-30-17 TIME 21:19:38

INTOXILYZER SERIAL # 80-001019 LOCATION Prescott YCSO

(x) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (x) 1. Clock time check.
- (x) 2. Date check.

OPERATIONAL TESTS

- (x) 1. Alcohol-free subject test result 0. 000 AC.
- (x) 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
- (x) 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
- (x) 4. Standard Calibration Check standard 0. 100 AC. Result 0. 097 AC.

Instrument is operating properly and accurately. YES X NO

COMMENTS

New Printer installed

SIGNATURE

[Signature]

Tank Serial # 27016100A7

Exp 11-5-2018

press 900 psi

DR 17-011168

INTOXILYZER 8000
Location: YCSO PRESCOTT
Serial Number: 80-001019
Core Version: 8105.48
03/30/2017 21:19:38

Standard Lot#: 27016100A1
Last Changed By: S. BLAKLEY SR #L6861

QAS: S. BLAKLEY SR #L6861
YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861
YAVAPAI CSO

Subject: SQAP, TEST
DOB: 10/20/1960
Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	21:20:58
Diagnostic Test	Pass	21:21:33
Air Blank	0.000	21:22:01
0.100 Cal Check	0.098	21:22:21
Air Blank	0.000	21:22:51
Subject Test	0.000	21:23:16
Air Blank	0.000	21:23:43
Wait	ABT*	
Air Blank	0.000	21:24:27

*Sequence Aborted

Not a Successfully Completed Test Sequence

DR 17-011168

INTOXILYZER 8000

Location: YCSO PRESCOTT

Serial Number: 80-001019

Core Version: 8105.48

03/30/2017

21:24:59

Standard Lot#: 27016100A1

Last Changed By: S. BLAKLEY SR #L6861

QAS: S. BLAKLEY SR #L6861

YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861

YAVAPAI CSO

Subject: SQAP, TEST

DOB: 10/20/1960

Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	21:26:08
Air Blank	0.000	21:26:37

*Improper Sample

Not a Successfully Completed Test Sequence

DR 17-011168

INTOXILYZER 8000
Location: YCSO PRESCOTT
Serial Number: 80-001019
Core Version: 8105.48
03/30/2017 21:26:55

Standard Lot#: 27016100A1
Last Changed By: S. BLAKLEY SR #L6861

QAS: S. BLAKLEY SR #L6861
YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861
YAVAPAI CSO

Subject: SQAP, TEST
DOB: 10/20/1960
Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	21:27:48
Air Blank	0.000	21:28:16

*RFI Detect

Not a Successfully Completed Test Sequence

DR 17-011168

INTOXILYZER 8000

Location: YCSO PRESCOTT

Serial Number: 80-001019

Core Version: 8105.48

03/30/2017

21:28:55

Standard Lot#: 27016100A1

Last Changed By: S. BLAKLEY SR #L6861

QAS: S. BLAKLEY SR #L6861
YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861
YAVAPAI CSO

Subject: SQAP, TEST

DOB: 10/20/1960

Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	21:30:04
Diagnostic Test	Pass	21:30:39
Air Blank	0.000	21:31:06
0.100 Cal Check	0.097	21:31:27
Air Blank	0.000	21:31:57
Subject Test	0.000*	21:35:04
Air Blank	0.000	21:35:36
Air Blank	0.000	21:36:03
Subject Test	XXX**	21:36:21
Air Blank	0.000	21:36:55

*Deficient Sample

**Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

DR 17-011168

COPY

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST S. Blakley Sr. AGENCY Yavapai County SO

DATE 3-30-17 TIME 21:40:26

INTOXILYZER SERIAL # 80-001019 LOCATION Prescott YCSO

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution NA AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0.098 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

DPS Form Exh G-2 (Rev 05-01)

Tank Serial # 27016100A1
Exp 11-5-2018
Press 900 Psi

DR17-011168

INTOXILYZER 8000

Location: YCSO PRESCOTT

Serial Number: 80-001019

Core Version: 8105.48

03/30/2017

21:40:26

Standard Lot#: 27016100A1

Last Changed By: S. BLAKLEY SR #L6861

QAS: S. BLAKLEY SR #L6861

YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861

YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	21:40:57
0.100 Cal Check	0.098	21:41:18
Air Blank	0.000	21:41:48

Preventative Maintenance
Performed -- Timer Reset