

COPY

EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000  
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. Murray AGENCY YCSO  
DATE 2/29/13 TIME 10:12:29  
INTOXILYZER SERIAL # 80-001455 LOCATION Budget 103

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution        AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (✓) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (✓) 4. Air blank completed.
- (✓) 5. Calibration check completed. Test results 0. 102 AC.
- (✓) 6. Air blank completed.
- (✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray NS6/5597

COPY

INTOXILYZER 8000

Location: YAVAPAI SO, BAGDAD

Serial Number: 80-001455

Core Version: 8105.48

07/29/2013

10:12:29

-----  
Standard Lot#: 13712100A2

Last Changed By: D. ASHBY #S62

-----  
QAS: J. MURRAY #N56

YAVAPAI YCSO

-----  
Operator: J. MURRAY #N56

YAVAPAI YCSO

-----  
Test                      g/210L                      Time  
-----  
Air Blank                      0.000                      10:13:26  
0.100 Cal Check                0.102                      10:13:47  
Air Blank                      0.000                      10:14:17

Preventative Maintenance  
Performed -- Timer Reset



EXHIBIT G-4  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST J. MURPHY AGENCY YCSO  
DATE 7/29/13 TIME 10:16:53  
INTOXILYZER SERIAL # 80-001455 LOCATION Breda 103

() 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- () 1. Clock time check.  
() 2. Date check.

OPERATIONAL TESTS

- () 1. Alcohol-free subject test result 0.000 AC.  
() 2. Error recognition logic system functioning.  
Not a Successfully Completed Test Sequence printed  
() 3. Proper sample recognition system.  
Not a Successfully Completed Test Sequence printed  
Deficient sample printed.  
() 4. Standard Calibration Check standard 0.100 AC. Result 0.101 AC.

Instrument is operating properly and accurately. YES  NO

COMMENTS 13712100A2 100 PSI  
7/1/2014

SIGNATURE J. Murphy NS6/5597



INTOXILYZER 8000

Location: YAVAPAI SO, BAGDAD

Serial Number: 80-001455

Core Version: 8105.48

07/29/2013

10:16:53

---

Standard Lot#: 13712100A2

Last Changed By: D. ASHBY

#S62

---

QAS: J. MURRAY

#N56

YAVAPAI YCSO

---

Operator: J. MURRAY

#N56

YAVAPAI CSO

---

Subject: SQAP, TEST

DOB: 10/20/1960

Sex: F

Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	10:18:08
Diagnostic Test	Pass	10:18:43
Air Blank	0.000	10:19:11
0.100 Cal Check	0.102	10:19:32
Air Blank	0.000	10:20:01
Subject Test	0.000	10:20:28
Air Blank	0.000	10:20:55
Wait	ABT*	
Air Blank	RFI**	10:22:14

\*Sequence Aborted

\*\*RFI Detect

**Not a Successfully Completed Test Sequence**



INTOXILYZER 8000  
Location: YAVAPAI SO, BAGDAD  
Serial Number: 80-001455  
Core Version: 8105.48  
07/29/2013

10:23:04

---

Standard Lot#: 13712100A2  
Last Changed By: D. ASHBY #S62

---

QAS: J. MURRAY #N56  
YAVAPAI YCSO

---

Operator: J. MURRAY #N56  
YAVAPAI CSO

---

Subject: SQAP, TES  
DOB: 10/20/1960  
Sex: F Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	IPS*	10:24:16
Air Blank	RFI**	10:24:23

\*Improper Sample  
\*\*RFI Detect

**Not a Successfully Completed Test Sequence**



INTOXILYZER 8000  
Location: YAVAPAI SO, BAGDAD  
Serial Number: 80-001455  
Core Version: 8105.48  
07/29/2013

10:25:32

---

Standard Lot#: 13712100A2  
Last Changed By: D. ASHBY #S62

---

QAS: J. MURRAY #N56  
YAVAPAI YCSO

---

Operator: J. MURRAY #N56  
YAVAPAI CSO

---

Subject: SQAP,TEST  
DOB: 10/20/1960  
Sex: F Weight: 115

---

15 Minute Deprivation Period? Yes

---

Test	g/210L	Time
Air Blank	0.000	10:26:27
Diagnostic Test	Pass	10:27:02
Air Blank	0.000	10:27:29
0.100 Cal Check	0.101	10:27:50
Air Blank	0.000	10:28:20
Subject Test	0.093*	10:31:27
Air Blank	0.000	10:32:00
Air Blank	0.000	10:32:27
Subject Test	XXX**	10:33:07
Air Blank	0.000	10:33:41

\*Deficient Sample

\*\*Invalid Sample

Begin new deprivation period

**Not a Successfully Completed Test Sequence**



EXHIBIT G-2  
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY  
STANDARD QUALITY ASSURANCE PROCEDURES  
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST J. Murray AGENCY YCSO  
DATE 7/29/13 TIME 10:34:04  
INTOXILYZER SERIAL # 80-001455 LOCATION Bagdad 103

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.  
OR  
Pour a standard alcohol concentration solution 2 AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach  $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0. 101 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE J. Murray NSB/5597

1. The first part of the document is a list of the names of the members of the committee.

2. The second part of the document is a list of the names of the members of the committee.

3. The third part of the document is a list of the names of the members of the committee.

INTOXILYZER 8000

Location: YAVAPAI SO, BAGDAD

Serial Number: 80-001455

Core Version: 8105.48

07/29/2013

10:34:04

-----  
Standard Lot#: 13712100A2

Last Changed By: D. ASHBY #S62  
-----

QAS: J. MURRAY #N56

YAVAPAI CSO  
-----

Operator: J. MURRAY #N56

YAVAPAI CSO  
-----

Test	g/210L	Time
Air Blank	0.000	10:34:36
0.100 Cal Check	0.101	10:34:57
Air Blank	0.000	10:35:27

Preventative Maintenance  
Performed -- Timer Reset

RECEIVED  
AUG 05 2013  
IC