

DR 18-034692

COPY

EXHIBIT G-2
 THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)
 ARIZONA DEPARTMENT OF PUBLIC SAFETY
 STANDARD QUALITY ASSURANCE PROCEDURES
 INTOXILYZER MODEL 8000
 STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Scott Blakley Sr. AGENCY Yavapai County SO
 DATE 9-19-18 TIME 22:45:27
 INTOXILYZER SERIAL # 80-01734 LOCATION YCSO Seligman

- (X) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
 OR
 Pour a standard alcohol concentration solution NA AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$.
- (X) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 4. Air blank completed.
- (X) 5. Calibration check completed. Test results 0.098 AC.
- (X) 6. Air blank completed.
- (X) 7. Remove instrument record. Attach the record to the completed checklist.

SIGNATURE Scott Blakley Sr.

DPS Form I-201 C (Rev. 03-01)

Tank Serial # 06717100A5

Exp 5-5-2019
 PRESS 900 PSI

Instrument Returned to service

DR 19-03-4692

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number: 80-001734
Core Version: 8105.48
09/19/2018 22:45:27

Standard Lot#: 0617100A5
Last Changed By: S. BLAKLEY SR #L6861

QAS: S. BLAKLEY SR #L6861
YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861
YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	22:45:58
0.100 Cal Check	0.098	22:46:18
Air Blank	0.000	22:46:48

Preventative Maintenance
Performed -- Timer Reset

DR18-034692

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST Scott Bialley Sr. AGENCY Yavapai County SO

DATE 9-19-18 TIME 22:47:56

INTOXILYZER SERIAL # 80-01734 LOCATION YCSO Seligman

(X) 1. Check, Read "PUSH BUTTON TO START"

DIAGNOSIS NOTES

- (X) 1. Clock time check.
- (X) 2. Date check.

OPERATIONAL TESTS

- (X) 1. Alcohol-free subject test result 0. 000 AC.
- (X) 2. Breath recognition logic system functioning.
Result: Successfully Completed Test Sequence printed
- (X) 3. Proper sample recognition system.
Result: Successfully Completed Test Sequence printed
No print sample printed.
- (X) 4. Standard Calibration Check standard 0. 100 AC. Result 0. 096 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

Instrument taken to service

SIGNATURE Scott Bialley Sr.

1 R: 039692

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number: 80-001734
Core Version: 8105.48
09/19/2018

22:47:56

Standard Lot#: 0617100A5
Last Changed By: S. BLAKLEY SR #L6861

CAS: S. BLAKLEY SR #L6861
YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861
YAVAPAI CSO

Subject: SOAP,TEST

DOB: 10/20/1960

Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	22:48:55
Diagnostic Test	Pass	22:49:30
Air Blank	0.000	22:49:58
0.100 Cal Check	0.096	22:50:19
Air Blank	0.000	22:50:48
Subject Test	0.000	22:51:04
Air Blank	0.000	22:51:32
wait	ABT*	
Air Blank	0.000	22:52:15

Sequence Aborted

Not a Successfully Completed Test Sequence

DR 18-034/692

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number: 80-001734
Core Version: 8105.48
09/19/2018

22:53:40

Standard Lot#: 0617100A5
Last Charged By: S. BLAKLEY SR #L6861

OAS: S. BLAKLEY SR #L6861
YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861
YAVAPAI CSO

Subject: SQAP, TEST
DOB: 10/20/1960

Sex: M Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	22:54:24
Air Blank	RFI**	22:54:36

Improper Sample
RFI Detect

Not a Successfully Completed Test Sequence

DR1803492

INTOXILYZER 8000

Location: YCSO SELIGMAN

Serial Number: 80-001734

Core Version: 8105.48

09/19/2018

22:55:24

Standard Lot#: 0617100A5

Last Changed By: S. BLAKLEY SR #L6861

OAS: S. BLAKLEY SR #L6861

YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861

YAVAPAI CSO

Subject: SQAP, TEST

DOB: 10/20/1960

Sex: F

Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
1st Blank	0.000	22:56:22
Diagnostic Test	Pass	22:56:57
2nd Blank	0.000	22:57:24
0.155 Cal Check	0.097	22:57:45
3rd Blank	0.000	22:58:15
Subject Test	0.000*	23:01:22
4th Blank	0.000	23:01:54
5th Blank	0.000	23:02:22
Subject Test	RGE**	23:02:45
6th Blank	0.000	23:03:23

Insufficient Sample

Range Exceeded

Not a Successfully Completed Test Sequence

Dr 8-034692

INFORMILYZER 8000
Location: YCSO SELIGMAN
Serial Number: 80-001734
Core Version: 8105.48
09/19/2018

23:04:17

Standard Lot#: 0617100A5
Last Changed By: S. BLAKLEY SR #L6861

QAS: S. BLAKLEY SR #L6861
YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861
YAVAPAI CSO

Subject: SQAP, TEST
DOB: 10/20/1960

Sex: F Weight: 115

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	23:05:14
Diagnostic Test	Pass	23:05:49
Air Blank	0.000	23:06:16
1000 Cal Check	0.097	23:06:37
Air Blank	0.000	23:07:07
Subject Test	0.295*	23:10:14
Air Blank	0.000	23:10:46
Air Blank	0.000	23:11:14
Subject Test	XXX**	23:11:33
Air Blank	0.000	23:12:06

Insufficient Sample

Invalid Sample

Begin new deprivation period

Not a Successfully Completed Test Sequence

DR 18-034692

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number: 80-001734
Core Version: 8105.48
09/19/2018

23:13:32

Standard Lot#: 0617100A5
Last Changed By: S. BLAKLEY SR #L6861

QAS: S. BLAKLEY SR #L6861
YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861
YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	23:14:01
0.100 Cal Check	0.096	23:14:22
Air Blank	0.000	23:14:52

Preventative Maintenance
Performed -- Timer Reset

DR 18-034692

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)
ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000
STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST Scott Blakely Sr. AGENCY Yavapai County SD
DATE 9-18-18 TIME 23:16:02
INTOXILYZER SERIAL # 90-01734 LOCATION YCSO Seligman

- (X) 1. Ensure the gas unit is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
2. Add a standard alcohol concentration solution NA AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.5^{\circ}\text{C}$.
- (X) 3. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- (X) 4. Go to the "Control Testing menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- (X) 5. Air blank completed.
- (X) 6. Calibration check completed. Test results 0.096 AC.
- (X) 7. All tests completed.

SIGNATURE [Signature]

INTOXILYZER 8000
Location: YCSO SELIGMAN
Serial Number: 80-001734
Core Version: 8105.48
09/19/2018

23:16:02

Standard Lot#: 0617100A5
Last Changed By: S. BLAKLEY SR #L6861

QAS: S. BLAKLEY SR #L6861
YAVAPAI CSO

Operator: S. BLAKLEY SR #L6861
YAVAPAI CSO

Test	g/210L	Time
Air Blank	0.000	23:16:31
0.100 Cal Check	0.096	23:16:52
Air Blank	0.000	23:17:22

Preventative Maintenance
Performed -- Timer Reset



ARIZONA DEPARTMENT OF PUBLIC SAFETY

INTOXILYZER MAINTENANCE REPORT

CR NO. I800017340006

Page 1 of 1

AGENCY: Yavapai County Sheriff's Office
Prescott, AZ 86301
OFFICER: BLAKLEY, #6861
DATE: September 11, 2018



INSTRUMENT

80010734

INSTRUMENT LOCATION

Yavapai County Sheriff's Office

MAINTENANCE REQUESTED

Preventative Maintenance

MAINTENANCE PERFORMED

Preventative Maintenance Performed:

Calibration, Verification

JOSEPH M SLOWINSKI, #6101, Electronics Technician
Central Regional Crime Laboratory
2102 W. Encanto Blvd., Phoenix, AZ 85009
602-223-2159