

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST T. ARIZWABERGSR #252 AGENCY PRESCOTT P.D.
DATE 092409 TIME 1243
INTOXILYZER SERIAL # 80-000337 LOCATION PRESCOTT P.D.

- (✓) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach 34°C ± 0.2°C
- (✓) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
(✓) 3. Go to the "Control Testing Menu". Select D for dry control test or "W" for wet control test. After selection is made press ENTER.
(✓) 4. Air blank completed.
(✓) 5. Calibration check completed. Test results 0. 099 AC.
(✓) 6. Air blank completed.
(✓) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE 

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

09/24/2009

12:43:37


Standard Lot#: 818201I

Last Changed By: G. KOESTER #212

QAS: T. ATTENBERGER #252

PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	12:44:04
Cal Check	0.099	12:44:23
Air Blank	0.000	12:44:52



Operator's Signature

EXHIBIT G-4
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD QUALITY ASSURANCE PROCEDURE

QA SPECIALIST T. AMEN AGENCY Prescott P.D.
DATE 09/24/09 TIME 1246
INTOXILYZER SERIAL # 80-000337 LOCATION Prescott P.D.

(✓) 1. Display Reads "PUSH BUTTON TO START"

DIAGNOSTIC TESTS

- (✓) 1. Clock time check.
(✓) 2. Date check.

OPERATIONAL TESTS

- (✓) 1. Alcohol-free subject test result 0. 000 AC.
(✓) 2. Error recognition logic system functioning.
Not a Successfully Completed Test Sequence printed
(✓) 3. Proper sample recognition system.
Not a Successfully Completed Test Sequence printed
Deficient sample printed.
(✓) 4. Standard Calibration Check standard 0. 100 AC. Result 0. 099 AC.

Instrument is operating properly and accurately. YES NO

COMMENTS

ALL FUNCTIONING CORRECTLY.

SIGNATURE T. Amen

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number : 80-000337
09/24/2009 12:46:16

Standard Lot#: 818201I
Last Changed By: G. KOESTER #212

QAS: T. ATTENBERGER #252
PRESCOTT PD

Operator: T. ATTENBERGER #252
PRESCOTT PD

Subject: SQAP, TEST
DOB: 01/01/1900
Sex: M Weight: 300

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:47:30
Diagnostic Test	Pass	12:47:58
Air Blank	0.000	12:48:25
0.100 Cal Check	0.099	12:48:44
Air Blank	0.000	12:49:13
Subject Test	0.000	12:49:36
Air Blank	0.000	12:50:04
Five Minute Wait	ABT*	
Air Blank	0.000	12:50:38

*Sequence Aborted

A handwritten signature in black ink, appearing to be 'T. Attenberger', with the number '14252' written to the right of the signature.

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number : 80-000337
09/24/2009 12:51:10

Standard Lot#: 818201I
Last Changed By: G. KOESTER #212

QAS: T. ATTENBERGER #252
PRESCOTT PD

Operator: T. ATTENBERGER #252
PRESCOTT PD

Subject: SQAP, TEST
DOB: 01/01/1900
Sex: M Weight: 300

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	IPS*	12:51:49
Air Blank	0.000	12:52:17

*Improper Sample

A handwritten signature in black ink, appearing to read "A. K. 2009" or similar, written in a cursive style.

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number : 80-000337
09/24/2009

12:52:57

Standard Lot#: 818201I
Last Changed By: G. KOESTER #212

QAS: T. ATTENBERGER #252
PRESCOTT PD

Operator: T. ATTENBERGER #252
PRESCOTT PD

Subject: SQAP, TEST
DOB: 01/01/1900
Sex: M Weight: 300

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	RFI*	12:53:45
Air Blank	0.000	12:54:14

*RFI Detect

A handwritten signature in black ink, appearing to read "O. A. [unclear]". The signature is written in a cursive style with a large initial "O" and a long horizontal stroke.

Not a Successfully Completed Test Sequence

INTOXILYZER 8000
Location: PRESCOTT P.D
Serial Number : 80-000337
09/24/2009

12:54:47

Standard Lot#: 818201I
Last Changed By: G. KOESTER #212

QAS: T. ATTENBERGER #252
PRESCOTT PD

Operator: T. ATTENBERGER #252
PRESCOTT PD

Subject: SQAP, TEST
DOB: 01/01/1900
Sex: M Weight: 300

15 Minute Deprivation Period? Yes

Test	g/210L	Time
Air Blank	0.000	12:55:49
Diagnostic Test	Pass	12:56:16
Air Blank	0.000	12:56:44
0.100 Cal Check	0.099	12:57:02
Air Blank	0.000	12:57:32
Subject Test	0.026*	13:00:40
Air Blank	0.000	13:01:13
Air Blank	0.000	13:01:40
Subject Test	XXX**	13:02:07
Air Blank	0.000	13:02:42

*Deficient Sample

**Invalid Sample

Begin new deprivation period



Not a Successfully Completed Test Sequence

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST T. ARSBERGER #252 AGENCY Prescott P.D.
DATE 092409 TIME 1304
INTOXILYZER SERIAL # 80-000337 LOCATION Prescott P.D.

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution .100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0. 099 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE [Signature]

INTOXILYZER 8000

Location: PRESCOTT P.D

Serial Number: 80-000337

09/24/2009

13:04:09

Standard Lot#: 818201I

Last Changed By: G. KOESTER #212

QAS: T. ATTENBERGER #252

PRESCOTT PD

Test	g/210L	Time
Air Blank	0.000	13:04:37
Cal Check	0.099	13:04:55
Air Blank	0.000	13:05:25



Operator's Signature