

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY

STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST SCARIM, ROGER G. AGENCY COTTONWOOD PD

DATE 3/26/13 TIME 0745

INTOXILYZER SERIAL # 80-000933 LOCATION COTTONWOOD PD

- MM*) 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 0.100 AC.
OR
Pour a standard alcohol concentration solution _____ AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- MM*
MM) 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- MM*
MM) 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- MM*
MM) 4. Air blank completed.
- MM*
MM) 5. Calibration check completed. Test results 0.097 AC.
- MM*
MM) 6. Air blank completed.
- MM*
MM) 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE *ofc MM M.* 592

INTOXILYZER 8000

Location: COTTONWOOD PD

Serial Number: 80-000933

Core Version: 8105.48

03/26/2013

07:48:36

Standard Lot#: 13711100A1

Last Changed By: R. SCARIM #592

QAS: R. SCARIM #592

COTTONWOOD PD

Operator: R. SCARIM #592

COTTONWOOD PD

Test	g/210L	Time
Air Blank	RFI*	07:48:53
Air Blank	0.000	07:49:21

*RFI Detect

13-003633

INTOXILYZER 8000
Location: COTTONWOOD PD
Serial Number: 80-000933
Core Version: 8105.48
03/26/2013 07:50:02

Standard Lot#: 13711100A1
Last Changed By: R. SCARIM #592

QAS: R. SCARIM #592
COTTONWOOD PD

Operator: R. SCARIM #592
COTTONWOOD PD

Test	g/210L	Time
Air Blank	0.000	07:50:32
0.100 Cal Check	0.097	07:50:53
Air Blank	0.000	07:51:22

Preventative Maintenance
Performed -- Timer Reset

13-003633