

EXHIBIT G-2
THIS REPORT PREPARED PURSUANT TO DUTY IMPOSED BY A.A.C. R13-10-104 (A)

ARIZONA DEPARTMENT OF PUBLIC SAFETY
STANDARD QUALITY ASSURANCE PROCEDURES
INTOXILYZER MODEL 8000

STANDARD CALIBRATION CHECK PROCEDURE

QA SPECIALIST W. Lopez ZO AGENCY SEONA PD.
DATE 1-10-13 TIME 1250
INTOXILYZER SERIAL # 80-000924 LOCATION SEONA P.D.

- () 1. Ensure that gas tank is attached to instrument and contains a standard alcohol concentration solution 100 AC.
OR
Pour a standard alcohol concentration solution N/A AC, into a clean dry simulator and assemble the simulator. Ensure that a tight seal has been made. Turn on the simulator and allow temperature to reach $34^{\circ}\text{C} \pm 0.2^{\circ}\text{C}$
- () 2. Intoxilyzer 8000 display reads "PUSH BUTTON TO START"
- () 3. Go to the "Control Testing Menu". Select "D" for dry control test or "W" for wet control test. After selection is made press ENTER.
- () 4. Air blank completed.
- () 5. Calibration check completed. Test results 0.099 AC.
- () 6. Air blank completed.
- () 7. Remove printed record. Attach the record to the completed checklist.

SIGNATURE W. Lopez ZO

DPS Form Exh G-2 (Rev 05-01)

PASS

I hereby certify that the above and forgoing is a true and correct copy of the record of periodic maintenance for the Intoxilyzer 8000 referenced therein maintained by the Seona Police Department pursuant to the requirements of the Arizona Department of Public Safety.

Date this 10TH day of JAN., 2013

W. Lopez ZO
Quality Assurance Specialist

INTOXILYZER 8000
Location: SEDONA
Serial Number: 80-000924
Core Version: 8105.48
01/10/2013

12:51:48

Standard Lot#: 09312100A1
Last Changed By: W. LOPEZ #30

QAS: W. LOPEZ #30
SEDONA PD

Operator: W. LOPEZ #30
SEDONA PD

Test	g/210L	Time
Air Blank	0.000	12:52:24
0.100 Cal Check	0.099	12:52:44
Air Blank	0.000	12:53:14

Preventative Maintenance
Performed -- Timer Reset

I hereby certify that the above and forgoing is a true and correct copy of the record of periodic maintenance for the intoxilyzer 8000 referenced therein maintained by the Sedona Police Department pursuant to the requirements of the Arizona Department of Public Safety.

Date this 10TH day of JAN. 2013

W. Lopez #30
Quality Assurance Specialist